



Recombinant *Saccharomyces cerevisiae* Eukaryotic translation initiation factor 2 subunit alpha (SUI2)

Product Code	CSB-YP323232SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P20459
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSTSHCRFYE NKYPEIDDIV MVNVQQAEM GAYVKLLEYD NIEGMILLSE LSRRRIRSIQ KLIRVGKNDV AVVLRVDKEK GYIDLSKRRV SSEDIKCEE KYQKSKTVHS ILRYCAEKFK IPEELYKTI AWPLSRKFGH AYEAFKLSII DETVWEGIEP PSKDVLDELK NYISKRLTPQ AVKIRADVEV SCFSYEGIDA IKDALKSAED MSTEQMQVKV KLVAAPLYVL TTQALDKQKG IEQLESAIEK ITEVITKYGG VCNITMPPKA VTATEDAELQ ALLESKELDN RSDSEDDDEDE SDDE
Source	Yeast
Target Names	SUI2
Protein Names	Recommended name: Eukaryotic translation initiation factor 2 subunit alpha Short name= eIF-2-alpha
Expression Region	1-304
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.