



Recombinant *Drosophila melanogaster* Glutamine synthetase 1, mitochondrial (Gs1)

Product Code	CSB-BP323235DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P20477
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	GNT TRSQFLANSP NTALDKSILQ RYRNLETPAN RVQATYLWID GTGENIRLKD RVLDPKVPSSV EDLPDWQYDG SSTYQAHGEN SDTTLKPRAI YRDPFKPGKN DVIVLCDTYS ADGKPTASNK RAAFQAAIDL ISDQEPWFGI EQEYTLLEDVD GRPFGWPENG FPAPQGPYYC GVGADRYYAR DLVEAHVVAC LYAGIDFAGT NAEVMPAQWE FQIGPAGIKA CDDLWVSRYY LQRIAEYGV VVTFDPKPME GQWNGAGAHT NFSTKEMRAD GGKAIIEEAI EKLSKRHERH IKAYDPKEGK DNERRLVGRL ETSSIDKFSW GVANRAVSVR VPRGVATAGK GYLEDRRPSS NCDPYAVCNA IVRTCLLNE
Source	Baculovirus
Target Names	Gs1
Protein Names	Recommended name: Glutamine synthetase 1, mitochondrial EC= 6.3.1.2 Alternative name(s): Glutamate--ammonia ligase 1
Expression Region	28-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.