



Recombinant Salmonella typhimurium Probable adenosine monophosphate-protein transferase fic (fic)

Product Code	CSB-EP323274SXB-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P20751
Product Type	Recombinant Protein
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Purity	>85% (SDS-PAGE)
Sequence	MSDKFGEGRD PYLYPGLNVM RNRLGIHQAAQ RLAQAAYEMT ALRAATIELG PLVRGLPHLC AIHRQLYQDI FDWAGQLREV DIYQQGDTRFC HFAYIEKEGN ALMQDLEEEG YLVGLAHEKF VERLAHYCE INVLHPFRFG SGLAQRIFFE QLALHAGYAL SWQGIAVETW KQANQSGAMG DLSALRAIFQ KAISEARETE
Source	E.coli
Target Names	fic
Protein Names	Recommended name: Probable adenosine monophosphate-protein transferase fic EC= 2.7.7.n1 Alternative name(s): Cell filamentation protein fic
Expression Region	1-200
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.