



Recombinant Maize streak virus genotype A Replication-associated protein (C1/C2)

Product Code	CSB-MP323330MBK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P14978
Product Type	Recombinant Protein
Immunogen Species	Maize streak virus genotype A (isolate Nigeria) (MSV)
Purity	>85% (SDS-PAGE)
Sequence	MASSSSNRQF SHRNANTFLT YPKCPENPEI ACQMIWELVV RWIPKYILCA REAHKDGSLH LHALLQTEKP VRISDSRFFD INGFHPNIQS AKSVNRVRDY ILKEPLAVFE RGTFFIPRKSP FLGKSDSEVK EKKPSKDEIM RDIISHATSK EEYLSMIQKE LPFDWSTKLQ YFEYSANKLF PEIQEEFTNP HPPSSPDLLC NESINDWLQP NIFQSSDERS RKQSLYIVGP TRTGKSTWAR SLGVHNYWQN NVDWSSYNED AIYNIVDDIP FKFCPCWKQL VGCQRDFIVN PKYGKKKKVQ KKSKPTILA NSDEDWMKEM TPGQLEYFEA NCIYIMSPG EKWYSPPELP PTEAVHSDRS
Source	Mammalian cell
Target Names	C1/C2
Protein Names	Recommended name: Replication-associated protein Short name= Rep EC= 2.7.7.- EC= 3.1.21.-
Expression Region	1-360
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.