



Recombinant Human Uncharacterized protein C17orf66 (C17orf66)

Product Code	CSB-BP003369HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2RTY3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAYEKSTDIS DVSRSMFLYP WLEYPDKTKE LRKAMAPVHL PLSCYQMPKE EFPPSPECWR QHPSKPNSVP YCYFKKPEIY THWHDLYDQR EEREAEKMLR KMRDDCRYIK EVHQTHIKMF HLPMSKLTIK SEMRSRPLEP TQDPLKWQRL RELTKSLESP REDEQFYAAQ ALGCLRISDK FVMEALQQVA QTGPEKVKYE AYRTLAILGC LNKHVIRALI KQLKEKNEGQ RMETLTGLRM ALNSWAAVSK DKRTQVGDEG KLVPVLQTLI KKSSEASLE AALCLGFLRP CSNMVQEFLL QCLCQGLKTQ RMKALRMLVK VMHVHSAPVI KAILDQLCSS SVLEDRFEAT QMLKTIGLEQ IQAQGLEELT FNLLRRKTHN EPFLAVRQAV AQTVEELKLK PTMMNLVEAQ LMNPDATARQ EAVISLGVLG IRSPQVFHLL LDLLDAENHQ AVKKSLETL ILCASIDPWI QNKLKNKVLV VYEAPKTNVK AEPTRFQKEP ENPEELTIQD FRLAKLNPLF IAKSITKVGQ KKTPAFPPCC SKPRKHRPQV IGPWQPRIKK QLRVLAIEIAK
Source	Baculovirus
Target Names	HEATR9
Protein Names	Recommended name: Uncharacterized protein C17orf66
Expression Region	1-570
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.