



Recombinant Human Uncharacterized protein C17orf66 (C17orf66)

| | |
|--------------------------|--|
| Product Code | CSB-EP003369HU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | A2RTY3 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MAYEKSTDIS DVSRSMFLYP WLEYPDKTKE LRKAMAPVHL PLSCYQMPKE EFPPSPECWR QHPSKPNSVP YCYFKKPEIY THWHDLYDQR EEREAEKMLR KMRDDCRYIK EVHQTHIKMF HLPMSKLTIK SEMRSRPLEP TQDPLKWQRL RELTKSLESP REDEQFYAAQ ALGCLRISDK FVMEALQQVA QTGPEKVKYE AYRTLAILGC LNKHVIRALI KQLKEKNEGQ RMETLTGLRM ALNSWAAVSK DKRTQVGDEG KLVVPLQTLI KKSSEASLE AALCLGFLRP CSNMVQEFLL QCLCQGLKTQ RMKALRMLVK VMHVHSAPVI KAILDQLCSS SVLEDRFEAT QMLKTIGLEQ IQAQQLEELT FNLLRRKTHN EPFLAVRQAV AQTVEELKLK PTMMNLVEAQ LMNPDATARQ EAVISLGVLG IRSPQVFHLL LDLLDAENHQ AVKKSLETL ILCASIDPWI QNKLKNKVLV VYEAPKTNVK AEPTRFQKEP ENPEELTIQD FRLAKLNPLF IAKSITKVGQ KKTPAFPPCC SKPRKHRPQV IGPWQPRIKK QLRVLAIEIAK |
| Source | E.coli |
| Target Names | HEATR9 |
| Protein Names | Recommended name: Uncharacterized protein C17orf66 |
| Expression Region | 1-570 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |