



Recombinant Human cytomegalovirus Uncharacterized protein UL36 (UL36)

Product Code	CSB-EP323461HWV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P16767
Product Type	Recombinant Protein
Immunogen Species	Human cytomegalovirus (strain AD169) (HHV-5) (HCMV)
Purity	≥85% (SDS-PAGE)
Sequence	MDDLRLTLMA YGCIAIRAGD FNLNDFLEQ ECGTRLHVAW PERCFIQLRS RSALGPFVVGK MGTVCSQGAY VCCQEYLHPF GFVEGPGFMR YQLIVLIGQR GGIYCYDDL R DCVYELAPTM KDFLRNGFRH RDHFHTMRDY QRPMVQYDDY WNAVMLYRGD VESLSAEVTK RGYASYTIDD PFDECPDTHF AFWTHNTEVM KFKETSFSV V RAGGSIQTME LMIRTVPRIT CYHQLLGALG HEVPERKEFL VRQYVLVDTF GVVYGYDPAM DAVYRLAEDV VMFTCVMGKK GHRNHRFSGR REAIVRLEKT PTCQHPKKTP DPMIMFDEDD DDELSLPRNV MTHEEAESRL YDAITENLMH CVKLVTTDSP LATHLWPQEL QALCDSPALS LCTDDVEGVR QKLRARTGSL HHFELSYRFH DEDPETYMGF LWDIPSCDRC VRRRRFKVCD VGRRHIIPGA ANGMPPLTPP HVYMNN
Source	E.coli
Target Names	UL36
Protein Names	Recommended name: Uncharacterized protein UL36
Expression Region	1-476
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.