



Recombinant Human cytomegalovirus Virion egress protein UL53 (UL53)

Product Code	CSB-MP323464HWV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P16794
Product Type	Recombinant Protein
Immunogen Species	Human cytomegalovirus (strain AD169) (HHV-5) (HCMV)
Purity	>85% (SDS-PAGE)
Sequence	MSSVSGVRTP RERRSALRSL LRKRRQRELA SKVASTVNGA TSANNHGPEP SPADARPLT LHDLDIFRE HPELELKYLN MMKMAITGKE SICLPFNFHS HRQHTCLDIS PYGNEQVSRI ACTSCEDNRI LPTASDAMVA FINQTSNIMK NRNFYYGFCK SSELLKLSTN QPPIFQIYYL LHAANHDIVP FMHAEDGRLH MHVIFENPDV HIPCDCITQM LTAAREDYSV TLNIVRDHVV ISVLCHAVSA SSVKIDVTIL QRKIDEMDIP NDVSESFERY KELIQELCQS SGNNLYEEAT SSYAIRSPLT ASPLHVVSTN GCGPSSSSQS TPPHLHPPSQ ATQPHHYSHH QSQSQQHHHR PQSPPPPLFL NSIRAP
Source	Mammalian cell
Target Names	NEC1
Protein Names	Recommended name: Virion egress protein UL53 Alternative name(s): Primary envelopment factor UL53 Protein HFRF2
Expression Region	1-376
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.