



Recombinant Type I restriction enzyme Ecoprri (prD)

Product Code	CSB-BP323499ENL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P17224
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	≥85% (SDS-PAGE)
Sequence	MVDCTKPIAE SNNFIILDKY NPDWKITESY QSEGLEREL IQDLVNQGYE YLPTLNNTKA MLANVREQLQ NLNNVEFLEA EWRRFVETWM DKPSDGVVEK ARKIHHDDYVH DFVFDDGRIQ NIYLLDRKNI LRNKVQVIKQ FEQAGTHANR YDVTILVNGL PLVQIELKKR GVAIREAFNQ IHRYSKESFN SENSLFKYLQ LFVISNGTDT RYFANTTKRD KNSFDFTMNW AKSDNTLIKD LKDFTATFFQ KHTLLNVLVN YSVFDSSQTL LVMRPYQIAA TERILWKIKS SFTSEELVKT GKRWWYLAHY RFW
Source	Baculovirus
Target Names	prD
Protein Names	Recommended name: Type I restriction enzyme Ecoprri EC= 3.1.21.3
Expression Region	1-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.