



Recombinant Enterobacteria phage fr Assembly protein (A)

Product Code	CSB-MP323664ECK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P15966
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage fr (Bacteriophage fr)
Purity	>85% (SDS-PAGE)
Sequence	MRKFIPTERM SKSHVSVRE YADGELEDNS LPLIYRSNWS PGQYTSTGPR TKEWHYPSSY SRGAIGIKAL DQGKYARLGT SWGREFEERA GYGMSIDARS CYSLFPVSNQ LTWIDVPTNV ANRATTEVLG KVTQGNFNLG VALAEARSTA SQLSTQTIAL IKAYTAARRG NWRQALRYLA LNENRKFNSK SVASRWLELQ FGWMPLLSDI QGAYEMLT KV HLKAFMPMRA VRQVGQNVSL SGRLTSPAAS YKSTCNISRR IVIWFYINDA RLAWLSSLGI LNPLGIVWEK VPFSFLVDWL LPVGNMLEGL TAPIGCSYQS GTVTDVISGE STITADDIYG WDTVRPATAK VQISAVHRGV QSVWPTTGVY VKSPFSMVHT LDALALFRQR LWK
Source	Mammalian cell
Target Names	A
Protein Names	Recommended name: Assembly protein Alternative name(s): Maturation protein Short name= A protein
Expression Region	1-393
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.