



Recombinant Bovine Myocyte-specific enhancer factor 2A (MEF2A)

| | |
|--------------------------|---|
| Product Code | CSB-EP013670BO-B |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | A2VDZ3 |
| Product Type | Recombinant Protein |
| Immunogen Species | Bos taurus (Bovine) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MGRKKIQITR IMDERNRQVT FTKRKFGMLK KAYELSVLCD CEIALIIFNS SNKLFQYAST DMDKVLLKYT EYNEPHESRT NSDIVEALNK KEHRGCDSPD PDTSYVLTPH TEEKYKKINE EFDNMMRNHK IAPGLPPQNF SMSVTVPVTS PSALSYTNPG SSLVSPSLAA SSALADTSM L SPPQATLHRN VSPGAPQRPP STGSAGGMLS TSDLTVPNGA GSSPVGNGFV NSRASP NLIG TTGANS LGKV MPTKSPPPPG GGSLGMNSRK PDLRVVIPPS SKGMMPLNT QRISSSQATQ PLATPVVSVT TPSLPPQGLV YSAMPTAYNT DYSLTSADLS ALQGFNSPGM LSLGQVSAWQ QHHLGQAALN SLVAGGQLSQ GSNLSINTNQ NINIKSEPI S PPRDRMTPSG FQQQQQPQP P P P P Q A P Q P Q P R Q E V G R S P V DSLSSSSSSY DGSDREDPRG DFHSPVVLGR PPNSDRESP SVKRM RMDAW VT |
| Source | E.coli |
| Target Names | MEF2A |
| Protein Names | Recommended name: Myocyte-specific enhancer factor 2A |
| Expression Region | 1-492 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |