



# Recombinant Shufflon-specific DNA recombinase (rci)

<b>Product Code</b>	CSB-EP323698ENL-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P16470
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPSPRIRKMS LSRALDKYLK TVSVHKKGHQ QEFYRSNVIK RYPIALRNMD EITTVDIATY RDVRLAEINP RTGKPIGTNT VRLELALLSS LFNIARVEWG TCRTNPVELV RPKKVSSGRD RRLTSSEERR LSRYFREKNL MLYVIFHLAL ETAMRQGEIL ALRWEHIDLR HGV AHLPETK NGH SRDVPLS RRARNFLQMM PVNLHGNVFD YTASGFKN AW RIATQRLRIE DLHFHDLRHE AISRFFELGS LNVMEIAAIS GHRSMNMLKR YTHLRAWQLV SKLDARRRQT QKVAAWFVPY PAHITTINEE NGQKAHRIE GDFDNLHVTA TTKEEAVHRA SEVLLRTLAI AAQKGERVPS PGALPVNDPD YIMICPLNPG STPL
<b>Source</b>	E.coli
<b>Target Names</b>	rci
<b>Protein Names</b>	Recommended name: Shufflon-specific DNA recombinase
<b>Expression Region</b>	1-384
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.