



Recombinant Replication protein (rep)

Product Code	CSB-MP323821CMB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P18016
Product Type	Recombinant Protein
Immunogen Species	Clostridium perfringens
Purity	≥85% (SDS-PAGE)
Sequence	MRNALKLNKK NYIDNIHSRS KGWITRSVID KKGYSQWHYK YAEKDLDM DENIYITLNT FYKPCRLEN IKELNTLFID LDYYKTGKTK DQVLMDLEKN YFNQSIPIPN YVIDSGRGM YLIWIINAVPS KALPLWKAVQ EYLYNQLKYF GADRQALDAT RILRVPGSIN SKSKTVVNIL DEYEYIDLR EIQNGFLPEL KPYERKKGRP SKINYIYRER SLYYGRIQDI IKLCELREYD LKGHRELIF LYRYYLCSFT EDIEKALNDV LELNSMFRQH LSEREVIRAT RSAERCYLDK NKQYKYKNET LIELLEITEE EQKYMTIIS KKEYKRRENI RGKKNYQEQL KAKGKATKKE ELNVLRRKKIK ALKEKGFKNK EITLMLEVPI KTLERHITYM KKNGLL
Source	Mammalian cell
Target Names	rep
Protein Names	Recommended name: Replication protein Alternative name(s): ORF5
Expression Region	1-406
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.