



Recombinant *Drosophila virilis* Homeobox protein rough (ro)

Product Code	CSB-EP323838DMP-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P18264
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila virilis</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MLQQNPSSPL HAEQFIKSSN STPTATPAPA SATATAPIHT VLATPQRPS PRQFFERLYG HLETRTHNSE SGDIDVGTHA ETPYQSDGGS ASSPDISISD ERVSLVSFPS YELHAPTPDY SAYGAIANAT LPPALSFFAF SADPHINAGF SAFLARRRRK EGRQRRQRTT FSTEQTLRLE VEFHRNEYIS RSRRFELAET LRLSETQIKI WFQNRRAKDK RIEKAQIDQH YRNFVWANGF MSTIMGQTSY ATSAPATSSA AATAAAMLGS SYAATPALL PSTLVNNARR ADVNYVDIEA QHQQQQHHH QVHPDQQQQQ QQQQRRSRF VPPPTLINS
Source	E.coli
Target Names	ro
Protein Names	Recommended name: Homeobox protein rough
Expression Region	1-340
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.