



Recombinant *Saccharomyces cerevisiae* Ty transcription activator TEC1 (TEC1)

Product Code	CSB-MP323847SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P18412
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	SLKEDDFGK DNSRNIESYT GRIFDVYIQK DSYSQSALDD MFPEAVVSTA ACVKNEAEDN INLIDTHPQF ELVNTGLGAK SDDLKSPSAK ATFTDKQRKN EVPNISVSNY FPGQSSETSS TTESWTIGCD KWSEKVEEAF LEALRLIMKN GTTKIKIRNA NFGRNELISL YIKHKTNEFR TTKQISSHIQ VWKKTIQNKI KDSLTLSSKE KELLHLIEHG AEQTTENSNL FYDIFEEIID SLPSVSDSGS LTPKNLYVSN NSSGLSVHSK LLTPITASNE KKIENFIKTN AASQAKTPLI YAKHIYENID GYKCVPSKRP LEQLSPTELH QGDRPNKASF SNKKAILESA KKIEIEQRKI INKYQRISRI QEHESNPEFS SNSNSGSEYE SEEEVVPRSA TVTQLQSRPV PYYKNNGMPY SLSKVRGRPM YPRPAEDAYN ANYIQGLPQY QTSYFSQLLL SSPQHYEHSP HQRNFTPSNQ SHGNFY
Source	Mammalian cell
Target Names	TEC1
Protein Names	Recommended name: Ty transcription activator TEC1
Expression Region	2-486
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.