



Recombinant Enterobacteria phage BF23 Protein A2-A3 (A2-A3)

Product Code	CSB-MP323925ECE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P19348
Product Type	Recombinant Protein
Immunogen Species	Escherichia phage Bf23 (Enterobacteria phage BF23)
Purity	>85% (SDS-PAGE)
Sequence	MANVKTEKSA KFAWNEENTQ QAVTMYQQLI NENGLFANS DGLKEIAKTV GAASPVSVPS KLTSAKAYQK SDKPRKVGSG SSIRKAHYVR VIAKHAIKSG IIKDADDLAS LESAKLETLD AVAQLLGVD EVKQAAGE
Source	Mammalian cell
Target Names	A2-A3
Protein Names	Recommended name: Protein A2-A3
Expression Region	1-138
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.