



Recombinant Enterobacteria phage T4 Head vertex protein Gp24 (24)

Product Code	CSB-MP323981EDZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P19896
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Purity	>85% (SDS-PAGE)
Sequence	MAKINELLRE STTTNSNSIG RPNLVALTRA TTKLIYSDIV ATQRTNQPVA AFYGIKYLNP DNEFTFKTGA TYAGEAGYVD REQITELTEE SKLTLNKGDL FKYNNIVYKV LEDTPFATIE ESDLELALQI AIVLLKVRLF SDAASTSKFE SSDSEIADAR FQINKWQTAV KSRKLKTGIT VELAQDLEAN GFDAPNFLED LLATEMADEI NKDILQSLIT VSKRYKVTGI TDSGFIDLSY ASAPEAGRSL YRMVCEMVSH IQKESTYTAT FCVASARAAA ILAASGWLKH KPEDDKYLSQ NAYGFLANGL PLYCDTNSPL DYVIVGVVEN IGEKEIVGSI FYAPYTEGLD LDDPEHVGAF KVVVDPESLQ PSIGLLVRYA LSANPYTVAK DEKEARIIDG GDMDKMAGRS DLSVLLGVKL PKIIIDE
Source	Mammalian cell
Target Names	24
Protein Names	Recommended name: Head vertex protein Gp24
Expression Region	1-427
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.