



Recombinant Mycobacterium kansasii Diacylglycerol acyltransferase/mycolyltransferase Ag85B (fbpB)

Product Code	CSB-MP324138MKZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P21160
Product Type	Recombinant Protein
Immunogen Species	Mycobacterium kansasii
Purity	>85% (SDS-PAGE)
Sequence	FSRPGLPVEY LQVPSAAMGR SIKVQFQSGG DNPAVYLLD GLRAQDDYNG WDINTPAFEW YYQSGLSVIM PVGGQSSFYS DWYSPACGKA GCTTYKWETF LTSELPQWLS ANRSVKPTGS AAVGISMAGS SALILSVYHP QQFIYAGSLS ALMDPSQGMG PSLIGLAMGD AGGYKASDMW GPSSDPAWQR NDPSLHIPEL VANNTRLWIY CGNGTPSELG GANVPAEFLE NFVRSSNLKF QDAYNAAGGH NAVFNLDANG THSWEYWGAQ LNAMKGDLQA SLGAR
Source	Mammalian cell
Target Names	fbpB
Protein Names	Recommended name: Antigen 85-B Alternative name(s): Antigen 85 complex B Short name= 85B Short name= Ag85B Extracellular alpha-antigen Fibronectin-binding protein B Mycolyl transferase 85B EC= 2.3.1.-
Expression Region	41-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.