



Recombinant Arabidopsis thaliana Protochlorophyllide reductase B, chloroplastic (PORB)

Product Code	CSB-MP324147DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P21218
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	QTAA TSSPTVTKSV DGKKT LRKGN VVVTGASSGL GLATAKALAE TGKWNVIMAC RDFLKAERAA KSVGMPKDSY TVMHLDLASL DSVRQFVDNF RRTETPLDVL VCNAAVYFPT AKEPTYSAEG FELSVATNHL GHFLLARLLL DDLKKS DYPS KRLIIVGSIT GNTNTLAGNV PPKANLGLDR GLAGGLNGLN SSAMIDGGDF DGAKAYKDSK VCNMLTMQEF HRRFHEETGV TFASLYPGCI ASTGLFREHI PLFRALFPPF QKYITKGYVS ETESGKRLAQ VVSDPSLTKS GVYWSWNNAS ASFENQLSEE ASDVEKARKV WEISEKLVGL A
Source	Mammalian cell
Target Names	PORB
Protein Names	Recommended name: Protochlorophyllide reductase B, chloroplastic Short name= PCR B EC= 1.3.1.33 Alternative name(s): NADPH-protochlorophyllide oxidoreductase B Short name= POR B
Expression Region	67-401
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.