



Recombinant HTH-type transcriptional regulator rafR (rafR)

Product Code	CSB-BP324276ENL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P21867
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	>85% (SDS-PAGE)
Sequence	SLKAIATTL GISVTTVSRA LGGFSDVAAS TRERVEAEAR RRGYRPNTQA RRLKTGKTDA IGLVYPENDV PFNSGVFMDM VSCISRELAY HDIDLLLIAD DEHADCHSYM RLVESRRIDA LIIAHTLDDD PRITHLHKAG IPFLALGRVP QGLPCAWFDF DNHAGTWQAT QKLIALGHKS IALLSENTSH SYVIARRQGW LDALHEHGLK DPLRLVSPT RRAGYLAVME LMSLPAPPTA IITDNDLSGD GAAMALQLRG RLSGKEAVSL VVYDGLPQDS IIELDVA AVI QSTRSLVGRQ ISDMVYQIIN GASPELQIT WTPIFYPGST VHSPSF
Source	Baculovirus
Target Names	rafR
Protein Names	Recommended name: HTH-type transcriptional regulator rafR Alternative name(s): Raffinose operon repressor
Expression Region	2-336
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.