



Recombinant Abutilon mosaic virus Capsid protein (AR1, AV1)

Product Code	CSB-EP324289AAE-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P21942
Product Type	Recombinant Protein
Immunogen Species	Abutilon mosaic virus (isolate West India) (AbMV)
Purity	>85% (SDS-PAGE)
Sequence	MPKRDLPWRS MPGTSKTSRN ANYSPRARIG PRVDKASEWV HRPMYRKPRI YRTLRTADVP RGCEGPCKVQ SYEQRHDISH VGKVMCISDV TRNGITHRV GKRFCVKSVY ILGKIWM DEN IKLQNHTNSV MFWLVRDRRP YGTPMDFGHV FNMFDNEPST ATVKNDLRDR YQVLHKFYGK VTGGQYASNE QAIVKRFWKV NNHVVYNHQE AGKYENHTEN ALLLYMACTH ASNPVYATLK IRIYFYDSL N
Source	E.coli
Target Names	AR1
Protein Names	Recommended name: Capsid protein Alternative name(s): Coat protein Short name= CP
Expression Region	1-251
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.