



# Recombinant *Saccharomyces cerevisiae* Isocitrate dehydrogenase [NADP], mitochondrial (IDP1)

<b>Product Code</b>	CSB-EP324291SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P21954
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	FSKI KVKQPVVELD GDEMTRIIWD KIKKKLILPY LDVDLKYYDL SVESRDATSD KITQDAAEAI KKYGVGIKCA TITPDEARVK EFNLHKMWKS PNGTIRNILG GTVFREPIVI PRIPRLVPRW EKPIIIGRHA HGDQYKATDT LIPGPGSLEL VYKPSDPTTA QPQTLKVYDY KGSGVAMAMY NTDESIEGFA HSSFKLAIK KLNLFSTKN TILKKYDGRF KDIFQEVYEA QYKSKFEQLG IHYEHLIDD MVAQMIKSKG GFIMALKNYD GDVQSDIVAQ GFGSLGLMTS ILVTPDGKTF ESEAAHGTVT RHYRKYQKGE ETSTNSIASI FAWSRGLLKR GELDNTPALC KFANILESAT LNTVQQDGIM TKDLALACGN NERSAYVTTE EFLDAVEKRL QKEIKSIE
<b>Source</b>	E.coli
<b>Target Names</b>	IDP1
<b>Protein Names</b>	Recommended name: Isocitrate dehydrogenase [NADP], mitochondrial Short name= IDH EC= 1.1.1.42 Alternative name(s): IDP NADP(+)-specific ICDH Oxalosuccinate decarboxylase
<b>Expression Region</b>	17-428
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.