



# Recombinant *Saccharomyces cerevisiae* Transcriptional repressor OPI1 (OPI1)

<b>Product Code</b>	CSB-MP324292SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P21957
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSENQRLGLS EEEVEAAEVL GVLKQSCRQK SQPSEDVSQA DKMPASESST TPLNILDRVS NKIISNVVTF YDEINTNKR LKSIGRLLDD DDDEHDDYDY NDDEFFTNKR QKLSRAIAKG KDNLKEYKLN MSIESKKRLV TCLHLLKLAN KQLSDKISCL QDLVEKEQVH PLHKQDGNAR TTTGAGEDET SSDEDDDDDEE FFDASEQVNA SEQSIVVKME VVGTVKKVYS LISKFTANSL PEPARSQVRE SLLNLPTNWF DSVHSTSLPH HASFYANCE EQKVEQQQQQ QQQQQQQQLL QQQLLQQQQQ KRNKDGDDSA SPSSVTANG KVLILAKESL EMVRNVMGVV DSTLGKAEW VKQKQEVKEM IRERFLQQQQ QYRQQQQKDG NYVKPSQDNV DSKD
<b>Source</b>	Mammalian cell
<b>Target Names</b>	OPI1
<b>Protein Names</b>	Recommended name: Transcriptional repressor OPI1 Alternative name(s): Negative regulator of phospholipid biosynthesis Overproducer of inositol protein 1
<b>Expression Region</b>	1-404
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.