



Recombinant Escherichia coli Protein mcrC (mcrC)

Product Code	CSB-MP324437ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P15006
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MEQPVIPVRN IYYMLTYAWG YLQEIKQANL EAIPGNLLD ILGYVLNKGV LQLSRRGLEL DYNPNTEIIP GIKGRIEFAK TIRGFHLNHG KTVSTFDMLN EDTLANRIIK STLAILIKHE KLNSTIRDEA RSLYRKLPGI STLHLTPQHF SYLNGGKNTR YYKFVISVCK FIVNNSIPGQ NKGHYRFYDF ERNEKEMSL YQKFLYEFGR RELTSANTTR SYLKWDASSI SDQSLNLLPR METDITIRSS EKILIVDAKY YKSIFSRRMG TEKFHSQONLY QLMNYLWSLK PENGENIGGL LIYPHVDTAV KHRYKINGFD IGLCTVNLGQ EWPCIHQELL DIFDEYLK
Source	Mammalian cell
Target Names	mcrC
Protein Names	Recommended name: Protein mcrC
Expression Region	1-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.