



Recombinant Escherichia coli Protein traH (traH)

Product Code	CSB-EP324441ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P15069
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	DVNSDM NQFFNKLGFASNTTQPGVWQ GQAAGYAYGG SLYARTQVKN VQLISMTLPD INAGCGGIDA YLGSFSFING EQLQRFVKQI MSNAAGYFFD LALQTTVPEI KTAKDFLQKM ASDINSMNLS SCQAAQGIIG GLFPRTQVVSQ QKVCQDIAGE SNIFADWAAS RQGCTVGGKS DSVRDKASDK DKERVTKNIN IMWNALSKNR MFDGNKELKE FVMTLTGSLV FGPNGEITPL SARTTDRSII RAMMEGGTAK ISHCNDSKDC LKVVADTPVT ISRDNALKSQ ITKLLASIQN KAVSDTPLDD KEKGFISSTT IPVFKYLVDP QMLGVSNSMI YQLTDYIGYD ILLQYIQELI QQARAMVATG NYDEAVIGHI NDNMNDATRQ IAAFQSQVQV QQDALLVVDR QMSYMRQQLS ARMLSRQNN YHFGGSTL
Source	E.coli
Target Names	traH
Protein Names	Recommended name: Protein traH
Expression Region	25-458
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.