



Recombinant Pig Myeloid differentiation primary response protein MyD88 (MYD88)

Product Code	CSB-BP015282PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2TF48
Product Type	Recombinant Protein
Immunogen Species	Sus scrofa (Pig)
Purity	>85% (SDS-PAGE)
Sequence	MAAGGSGAES APPTPSMSSL PLAALNVRVR HRLSLFLNVR TQVAADWTGL AEEMNFEYLE IRRLETHPDP TRSLDDWQG RPGASVGRLL ELLAKLGRDD VLVELGPSIE EDCRKYILKQ QQEAAEKPLQ VDSVDSSIPW MSGITIRDDP LGQMPEHFDA FICYCPSDIQ FVQEMIRQLE QTNYRLKLCV SDRDVLPGTC VWSIASIELIE KRCRRMVVVV SDDYLQSKEC DFQTKFALSL SPGAHQKRLI PVKYKSMKKE FPSILRFITV CDYTNPCTKS WFWTRLARAL SLP
Source	Baculovirus
Target Names	MYD88
Protein Names	Recommended name: Myeloid differentiation primary response protein MyD88
Expression Region	1-293
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.