



Recombinant Human cytomegalovirus Uncharacterized protein UL116 (UL116)

Product Code	CSB-EP324575HWV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P16833
Product Type	Recombinant Protein
Immunogen Species	Human cytomegalovirus (strain AD169) (HHV-5) (HCMV)
Purity	≥85% (SDS-PAGE)
Sequence	VETNAT TVTSTTAAAA TTNTTVATTG TTTTSPNVTS TTSNTVTTPT TVSSVSNLTS STTSIPISTS TVSGTRNTGN NNTTTIGTNA TSPSPSVSIL TTVTPAATST ISVDGVVTAS DYTPTFDDLE NITTTTRAPTR PPAQDLCSHN LSIILYEEES QSSVDIAVDE EEPELEDDDE YDELWFLPYF EAECNRNYTL HVNHS CDYSV RQSSVSFPPW RDIDSVTFVP RNLSNCSAHG LAVIVAGNQT WYVNPFLSLAH LLDAIYNVLG IEDLSANFRR QLAPYRHTLI VPQT
Source	E.coli
Target Names	UL116
Protein Names	Recommended name: Uncharacterized protein UL116
Expression Region	25-314
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.