



Recombinant Human cytomegalovirus 65 kDa phosphoprotein (UL83), partial

Product Code	CSB-EP324690HWX-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P18139
Product Type	Recombinant Protein
Immunogen Species	Human cytomegalovirus (strain Towne) (HHV-5) (Human herpesvirus 5)
Purity	>85% (SDS-PAGE)
Source	E.coli
Target Names	UL83
Protein Names	Recommended name: 65 kDa phosphoprotein Short name= pp65 Alternative name(s): 64 kDa matrix phosphoprotein Short name= pp64 GP64 Phosphoprotein UL83 Tegument protein UL83
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.