



Recombinant *Pectobacterium carotovorum* subsp. *carotovorum* Endo-polygalacturonase (peh)

Product Code	CSB-EP324694ESP-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P18192
Product Type	Recombinant Protein
Immunogen Species	<i>Pectobacterium carotovorum</i> subsp. <i>carotovorum</i> (<i>Erwinia carotovora</i> subsp. <i>carotovora</i>)
Purity	>85% (SDS-PAGE)
Sequence	AFASDSR TVSEPKAPSS CTVLKADSST ATSTIQKALN NCGQGKAVKL SAGSSSVFLS GPLSLPSGVS LLIDKGVTLR AVNNAKSFEN APSSCGVVDT NGKGCDAFIT ATSTTNSGIY GPGTIDGQGG VKLQDKKVSF WDLAADAQVK KLKQNTPRLI QINKSKNFTL YNVSLINSPN FHVVFSDGDG FTAWKTTIKT PSTARNTDGI DPMSSKNITI AHSNISTGDD NVAIKAYKGR SETRNIILH NEFGTGHGMS IGSETMGVYN VTVDDLIMTG TTNGLRIKSD KSAAGVVNGV RYSNVVMKNV AKPIVIDTVY EKKEGSNVPD WSDITFKDIT SQTKGVVVLN GENAKKPIEV TMKNVCLTSD STWQIKNVTV KK
Source	E.coli
Target Names	peh
Protein Names	Recommended name: Endo-polygalacturonase EC= 3.2.1.15
Expression Region	24-402
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.