



Recombinant *Saccharomyces cerevisiae* Regulator of drug sensitivity 2 (RDS2)

Product Code	CSB-MP324820SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P19541
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSANSGVKRA SKAFKTCLFC KRSHVVC DKQ RPCSRCVKRD IAHLCREDDI AVPNEMPSQH ESSPNDNNIQ GKYANKAHTG IPSDYQNEPV NKSGSTYGEE LSPKLDSSLV NDTTSLLLPQ QPVFVSENVG SEFSSLNEFL SMLNPLLTQ TSLSSSSASN VHLENGSQTT QSPLEYQNDN RRDEIGVARQ ENRSPTIMSG SSNSISKGDK QDQEKEESRI LANANENSAP TPKEQFFLTA ADPSTEMTPE HRLKLVINAK LEAGLLKPYN YAKGYARLQD YMDKYMNQSS KQRILKPLST IRPAFRTIAR SLKDVDLVLV EESFERMLLS YDRVFTSM SM PACLRRTGE IYRANKEFAS LVDCTVDDL R DGKLAIYELM TEESAVNFW E KYGSI AFDKG QKAVLTSCSL RTKDGIRKRP CCFSTIRRD RYNIPICIVG NFIPLS
Source	Mammalian cell
Target Names	RDS2
Protein Names	Recommended name: Regulator of drug sensitivity 2
Expression Region	1-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.