



Recombinant Bacillus subtilis Homoserine dehydrogenase (hom)

Product Code	CSB-BP324827BRJ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P19582
Product Type	Recombinant Protein
Immunogen Species	Bacillus subtilis (strain 168)
Purity	>85% (SDS-PAGE)
Sequence	MKAIRVGLLG LGTVGSGVVK IIQDHQDKLM HQVGCPVTIK KVLVKDLEKK REVDLPKEVL TTEVYDVIDD PDVDVVIEVI GGVEQTKQYL VDALRSKKHV VTANKDLMAV YGSELLAEAK ENGCDIYFEA SVAGGIPILR TLEEGLSSDR ITKMMGIVNG TTNFILTKMI KEKSPYEEVL KEAQDLGFAE ADPTSDVEGL DAARKMAILA RLGFSMNVDL EDVKVKGISQ ITDEDISFSK RLGTYMKLIG IAQRDGSKIE VSVQPTLLPD HHPLSAVHNE FNAVYVYGEA VGETMFYGPG AGSMPTATSV VSDLVAVMKN MRLGVTGNSF VGPQYEKNMK SPSDIYAQQF LRIHVKDEVG SFSKITSVFS ERGVSEKIL QLPIKGHDEL AEIVIVTHHT SEADFSDILQ NLNDLEVVQE VKSTYRVEGN GWS
Source	Baculovirus
Target Names	hom
Protein Names	Recommended name: Homoserine dehydrogenase Short name= HDH EC= 1.1.1.3
Expression Region	1-433
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.