



Recombinant Human PRAME family member 22 (PRAMEF22)

| | |
|--------------------------|---|
| Product Code | CSB-EP018618HU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | A3QJZ6 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MRMSLQAPRR LLELAGQSLI GDQALAISIL DELPRELFPP LFVEAFTSRR CEVLKVMVQA WPFPCPLPGS LMKTPDLEIL HYVVDGIDCL LAQKVRPRRW KLQVLELRDV DENFWTIWSG ARPLSCSPEA MSKRQTVEDC PRTGEKQPLK VFMDVCLKEK FMDEDLSFFS GWVQHRRGSV HLCCTKVVNY SMSILNFRNI LETVYPDSIQ VLEIWNMCWP CMIVEFSRYL SQMRNLRKLF ISDGCYLLS SDSQEQLVAE FSSVLLRLEY LQMLYVRRVC FFRGHLDQLI RCLRSPLETL ALTYGFLEKV DLKCLPRYPS LSQKQLNLS HGALRFIRLE PLRALLEKVA ATLQTLFLVD CGIRDSKLRV ILPALSCSN LTTFCFHGND TSMDGLKDLL RHTGRLSNLS LETYPAPRES LDDRGRVISE LLTPLQAELM RILREVREP RIFFGPVSCP CCGTSPTEQL EFNFLWGRP A |
| Source | E.coli |
| Target Names | PRAMEF22 |
| Protein Names | Recommended name: PRAME family member 22 |
| Expression Region | 1-481 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |