



Recombinant Human Protein LCHN (LCHN)

Product Code	CSB-MP012231HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A4D1U4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVEQGDAAPL LRWAEGPAVS LPQAPQPQAG GWGRGGGGGA RPAAEPPRRR EPEEPAAPEV LLQPGRLELG DVEEDQVVAV FVVTFDPRSG NMVEWCLPQD IDLEGVEFKS MASGSHKIQS DFIYFRKGP FGLACFANMP VESELERGAR MKSVGILSPS YTLRYMHF LENQVRHQLE MPGHYSHLAA FYEDKKGVLH AGPGRGSSLP PVYWLPSIHR YMYPEMKITH PAGCMSQFIK FFGEQILILW KFALLRKRIL IFSPPPVGVV CYRVYCCCCL ANVSLPGIGG TIPESKPFY VNVADIESLE VEVSYVACTT EKIFEEKREL YDVYVDNQNV KTHHDHLQPL LKINSADREK YRRLNEQRQM LLYSQEVEED YNPCEEDLFV LFFLEQNNRI FQTLLEVSAS QDKTLTAEHA RGMGLDPQGD RSFLDLLEA YGIDVMLVID NPCCP
Source	Mammalian cell
Target Names	LCHN
Protein Names	Recommended name: Protein LCHN
Expression Region	1-455
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.