



Recombinant Pseudomonas aeruginosa High-affinity branched-chain amino acid transport ATP-binding protein BraF (braF)

Product Code	CSB-MP325121EZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P21629
Product Type	Recombinant Protein
Immunogen Species	Pseudomonas aeruginosa (strain ATCC 15692 / DSM 22644 / CIP 104116 / JCM 14847 / LMG 12228 / 1C / PRS 101 / PAO1)
Purity	≥85% (SDS-PAGE)
Sequence	MSRPILEVSG LTMRFGGLLA VNGVNLKVEE KQVVSMIGPN GAGKTTVFNC LTGFYQPTGG LIRLDGEEIQ GLPGHKIARK GVVRTFQNVN LFKEMTAVEN LLVAQHRHLN TNFLAGLFKT PAFRRSERA MEYAAHWLEE VNLTEFANRS AGTLAYGQQR RLEIARCM MT RPRILMLDEP AAGLNPKETD DLKALIAKLR SEHNVTVLLI EHDMKLVMSI SDHIVVINQG APLADGTPEQ IRDNPDVIKA YLGEA
Source	Mammalian cell
Target Names	braF
Protein Names	Recommended name: High-affinity branched-chain amino acid transport ATP-binding protein BraF
Expression Region	1-255
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.