



Recombinant *Drosophila melanogaster* Chromosomal protein D1 (D1)

Product Code	CSB-EP325206DLU-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P22058
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MEEVAVKKRG RPSKASVGGK SSTAAVA AIS PGIKKRGRPA KNKGSSGGGG QRGRPPKASK IQNDEDPEDE GEEDGDGDGS GAELANNSSP SPTKGRGRPK SSGGAGSGSG DSVKTPGSAK KRKAGRPKKH QPSDSENE DD QDEDDDG NSS IEERRPVGRP SAGSVNLNIS RTGRGLGRPK KRAVESNGDG EPQVPKGRGR PPQNKSGSGG STGYVPTGRP RGRPKANAAP VEKHEDNDDD QDDENS GEEE HSSPEKTVVA PKKRGRPSLA AGKVSKEETT KPRSRPAKNI DDDADDADSA DQGQHNSKKE SNDEDRAVDG TPTKGDGLKW NSDGENDAND GYVSDNYNDS ESVA A
Source	E.coli
Target Names	D1
Protein Names	Recommended name: Chromosomal protein D1
Expression Region	1-355
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.