



Recombinant Candida albicans ADP-ribosylation factor (ARF1)

Product Code	CSB-YP325259CZD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P22274
Product Type	Recombinant Protein
Immunogen Species	Candida albicans (strain SC5314 / ATCC MYA-2876) (Yeast)
Purity	>85% (SDS-PAGE)
Sequence	GLTISKLFA SLLGRREMRI LMVGLDAAGK TTILYKCLKLG EIVTTIPTIG FNVETVEYKN ISFTVWDVGG QDKIRPLWRY YFQNTQGIIF VVDSNDRDRI NEAREELQSM LNEDELKDAV LLVLANKQDL PNAMNAAEIT EKMGLHSIMN RPWFIQATCA TTGDGLYEGL EWLSNQVGK
Source	Yeast
Target Names	ARF1
Protein Names	Recommended name: ADP-ribosylation factor
Expression Region	2-179
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.