



Recombinant Flavobacterium sp. N-acylmannosamine 1-dehydrogenase

Product Code	CSB-MP325300FBT
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P22441
Product Type	Recombinant Protein
Immunogen Species	Flavobacterium sp. (strain 141-8)
Purity	>85% (SDS-PAGE)
Sequence	TTAGVSRRP GRLAGKAAIV TGAAGGIGRA TVEAYLREGA SVVAMDLAPR LAATRYEEPG AIPIACDLAD RAAIDAAMAD AVARLGGLDI LVAGGALKGG TGNFLDLSA DWDRYVDVNM TGTFILTCRAG ARAMVAAGAG KDGRSARIIT IGSVNSFMAE PEAAAYVAAK GGVAMLTRAM AVDLARHGIL VNMIAPGPVD VTGNNTGYSE PRLAEQVLDE VALGRPGLPE EVATAAVFLA EDGSSFITGS TITIDGGLSA MIFGGMREGR R
Source	Mammalian cell
Protein Names	Recommended name: N-acylmannosamine 1-dehydrogenase Short name= NAM-DH EC= 1.1.1.233
Expression Region	2-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.