



Recombinant Solanum lycopersicum Suberization-associated anionic peroxidase 1 (TAP1)

Product Code	CSB-MP325320SZY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P15003
Product Type	Recombinant Protein
Immunogen Species	Solanum lycopersicum (Tomato) (<i>Lycopersicon esculentum</i>)
Purity	>85% (SDS-PAGE)
Sequence	RNTYE AIIMKNGSLL KNVSPDFDSL ESGVASILTL NNNKKRNSDK YLRQQLTPEA CVFSAVRVV DSAIDAETRM GASLIRLHFH DCFVDGCDGG ILLDDINGTF TGEQNPPNA NSARGYEVIA QAKQSVINTC PNVSVSCADI LAAARDSVA KLGGQTYVA LGRSDARTAN FSGAINQLPA PFDNLTVQIQ KFSDKNFTLR EMVALAGAHT VGFARCSTVC TSGNVNPAAQ LQCNC SATLT DSDLQQLDTT PTMFDKVVYD NLNSNQGIMF SDQVLTGDAT TAGFVTDYSN DVNVFLGDFA AAMIKMGDLP PSAGAQLAIR DVCSRVNPTS VASM
Source	Mammalian cell
Target Names	TAP1
Protein Names	Recommended name: Suberization-associated anionic peroxidase 1 EC=1.11.1.7 Alternative name(s): TMP1
Expression Region	26-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.