



Recombinant Cucumber necrosis virus RNA silencing suppressor p19 (ORF4)

Product Code	CSB-YP325335CUF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P15184
Product Type	Recombinant Protein
Immunogen Species	Cucumber necrosis virus (CNV)
Purity	>85% (SDS-PAGE)
Sequence	MERAIQRSDA REQANSERWD GRCGGTITPF KLPDESPSLLEWRLHNSEES EDKDHPLGFK ESWDFGKVVFKRYLRYDGTE TSLHRTLGSWERNSVNDAAAS RFLGVSQIGC TYSIRFRGSC LTLSGGSRTL QRLIEMAIRT KRTMLQLTPC EVEGNVSRRR PQGSEAFENK ESE
Source	Yeast
Target Names	ORF4
Protein Names	Recommended name: RNA silencing suppressor p19 Alternative name(s): 19 kDa symptom severity modulator
Expression Region	1-173
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.