



Recombinant Zea mays 30S ribosomal protein S2, chloroplastic (rps2)

Product Code	CSB-EP325398ZAX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P16037
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	>85% (SDS-PAGE)
Sequence	MTRRYWNINL KEMIEAGVHF GHGIKKWNP K MAPYISAKRK GTHITNLART ARFLSEACDL VFDAASQGKS FLIVGTTKRA ADLVASAAIR SRCHYVNKKW FSGMLTNWSI TKTRLSQFRD LRAEEKMGKF HHLPKRDAAI LKRKLSTLQR YLGGIKYMTR LPDIVIVLDQ QKEYIALQEC AILGIPTISL VDTNCDPDLA NISIPANDDT MTSIRLILNK LVFAISEGRS LYIRNR
Source	E.coli
Target Names	rps2
Protein Names	Recommended name: 30S ribosomal protein S2, chloroplastic
Expression Region	1-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.