



# Recombinant Escherichia coli L-serine dehydratase 1 (sdaA)

<b>Product Code</b>	CSB-EP325402ENV-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P16095
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MISLFDMPFKV GIGPSSSHTV GPMKAGKQFV DDLVEKGLLD SVTRVAVDVY GSLSLTGKGH HTDIAIIMGL AGNEPATVDI DSIPGFIRDV EERERLLLAQ GRHEVDFPRD NGMRFHNGNL PLHENGMIQH AYNGDEVVYS KTYYSIGGGF IVDEEHFGQD AANEVSVYPY FKSATELLAY CNETGYSLSG LAMQNELALH SKKEIDEYFA HVWQTMQACI DRGMNTEGVL PGPLRVPRRA SALRRMLVSS DKLSNDPMNV IDWVNMFALA VNEENAAGGR VVTAPTNGAC GIVPAVLAYY DHFIESVSPD IYTRYFMAAG AIGALYKMNA SISGAEVGCQ GEVGVAC SMA AAGLAELLGG SPEQVCVAAE IGMEHNLGLT CDPVAGQVQV PCIERNAIAS VKAINAARMA LRRTSAPRVS LDKVIETMYE TGKDMNAKYR ETSRGG LAIK VQCD
<b>Source</b>	E.coli
<b>Target Names</b>	sdaA
<b>Protein Names</b>	Recommended name: L-serine dehydratase 1 Short name= SDH 1 EC= 4.3.1.17 Alternative name(s): L-serine deaminase 1 Short name= L-SD1
<b>Expression Region</b>	1-454
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.