



Recombinant Drosophila melanogaster Guanine nucleotide-binding protein G (o) subunit alpha 47A

Product Code	CSB-YP325422DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P16378
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	GCAQSAEER AAAARSLIE RNLKEDGIQA AKDIKLLLLG AGESGKSTIV KQMIIHESG FTAEDFKQYR PVVYSNTIQS LVAILRAMPT LSIQYSNNER ESDAKMVFDV CQRMHDTEPF SEELLAAMKR LWQDAGVQEC FSRSNEYQLN DSAKYFLDDL DRLGAKDYQP TEQDILRTRV KTTGIVEVHF SFKNLNFKLF DVGGQRSERK KWIHCFEDVT AIFCVAMSE YDQVLHEDET TNRMQESLKL FDSICNNKWF TDTSIILFLN KKDLFEEKIR KSPLTICFPE YTGGQEYGEA AAYIQAQFEA KNKSTSKEIY CHMTCATDTN NIQVFVDAVT DVIIANNLRG CGLY
Source	Yeast
Target Names	Galphao
Protein Names	Recommended name: Guanine nucleotide-binding protein G(o) subunit alpha 47A
Expression Region	2-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.