



# Recombinant *Saccharomyces cerevisiae* Pyruvate dehydrogenase E1 component subunit alpha, mitochondrial (PDA1)

<b>Product Code</b>	CSB-YP325423SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P16387
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	ATLKTTD KKAPEDIEGS DTVQIELPES SFESYMLEPP DLSYETSKAT LLQMYKDMVI IRRMEMACDA LYKAKKIRGF CHLSVGQEI AVGIENAITK LDSIITSYRC HGFTFMRGAS VKAVLAELMG RRAGVSYGKG GSMHLYAPGF YGGNGIVGAQ VPLGAGLAFQA HQYKNEACD FTLYGDGASN QGQVFESFNM AKLWNLPPVF CCENNKYGMG TAASRSSAMT EYFKRGQYIP GLKVNMDIL AVYQASKFAK DWCLSGKGPL VLEYETYRYG GHMSMDPGTT YRTRDEIQHM RSKNDPIAGL KMHLIDLGIA TEAEVKAYDK SARKYVDEQV ELADAAPPPE AKLSILFEDV YVKGTTETPTL RGRIPEDTWD FKKQGFASRD
<b>Source</b>	Yeast
<b>Target Names</b>	PDA1
<b>Protein Names</b>	Recommended name: Pyruvate dehydrogenase E1 component subunit alpha, mitochondrial EC= 1.2.4.1 Alternative name(s): Pyruvate dehydrogenase complex component E1 alpha Short name= PDHE1-A
<b>Expression Region</b>	34-420
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.