



Recombinant Bovine F-actin-capping protein subunit alpha-1 (CAPZA1)

Product Code	CSB-EP004510BO-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A4FUA8
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	<p>ADFEDRVSD EEKVRIAAKF ITHAPPGEFN EVFNDVRLLL NNDNLLREGA AHAFAYQNM D QFTPVKIEGY EDQVLITEHG DLGNSRFLDP RNKISFKFDH LRKEASDPQP EEADGGLKSW RESCDSALRA YVKDHYSNGF CTVYAKTIDG QQTIIACIES HQFQPKNFWN GRWRSEWKFT ITPPTAQVVG VLKIQVHYE DGNVQLVSHK DVQDSVTVSN EAQTAKEFIK IIEHAENEYQ TAISENYQTM SDTTFKALRR QLPVTRTKVD WNKILSYKIG KEMQNA</p>
Source	E.coli
Target Names	CAPZA1
Protein Names	Recommended name: F-actin-capping protein subunit alpha-1 Alternative name(s): CapZ alpha-1
Expression Region	2-286
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	CAPZA1 is a member of the F-actin capping protein alpha subunit family. This gene encodes the alpha subunit of the barbed-end actin binding protein. The protein regulates growth of the actin filament by capping the barbed end of growing actin filaments.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.