



Recombinant Bovine F-actin-capping protein subunit alpha-1 (CAPZA1)

| | |
|--------------------------|---|
| Product Code | CSB-MP004510BO |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | A4FUA8 |
| Product Type | Recombinant Protein |
| Immunogen Species | Bos taurus (Bovine) |
| Purity | >85% (SDS-PAGE) |
| Sequence | ADFEDRVSD EEKVRIAAKF ITHAPPGEFN EVFNDVRLLL NNDNLLREGA AHAFAQYNMD QFTPVKIEGY EDQVLITEHG DLGNSRFLDP RNKISFKFDH LRKEASDPQP EEADGGLKSW RESCDSALRA YVKDHYSNGF CTVYAKTIDG QQTIACIES HQFQPKNFWN GRWRSEWKFT ITPPTAQVVG VLKIQVHYE DGNVQLVSHK DVQDSVTVSN EAQTAKEFIK IIEHAENEYQ TAISENYQTM SDTTFKALRR QLPVTRTKVD WNKILSYKIG KEMQNA |
| Source | Mammalian cell |
| Target Names | CAPZA1 |
| Protein Names | Recommended name: F-actin-capping protein subunit alpha-1 Alternative name(s): CapZ alpha-1 |
| Expression Region | 2-286 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Target Details | CAPZA1 is a member of the F-actin capping protein alpha subunit family. This gene encodes the alpha subunit of the barbed-end actin binding protein. The protein regulates growth of the actin filament by capping the barbed end of growing actin filaments. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |