



# Recombinant Enterobacteria phage T4 Baseplate structural protein Gp27 (27)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-BP325483EDZ   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | P17172  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Enterobacteria phage T4 (Bacteriophage T4)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MSMLQRPGYP NLSVKLFDSY DAWSNRNFVE LAATITTLTM RDSLYGRNEG<br>MLQFYDSKNI HTKMDGNEII QISVANANDI NNVKTRIYGC KHFSVSVDSK<br>GDNIIAIELG TIHSIENLKF GRPFFPDAGE SIKEMLGVIY QDRLLTPAI<br>NAINAYVPDI PWTSTFENYL SYVREVALAV GSDKFVFWQ DIMGVNMMDY<br>DMMINQEPYP MIVGEP SLIG QFIQELKYPL AYDFVWLTKS NPHKRDPMKN<br>ATYAHSFLD SSIPMITTGK GENSIVVRSR GAYSEMTYRN GYEEAIRLQT<br>MAQYDGYAKC STIGNFNLTG VVKIIFNSDK NQFKTEFYVD EVIHLSNNN<br>SVTHLYMFTN ATKLETIDPV KVKNEFKSDT TTESSSSNK Q |
| <b>Source</b>            | Baculovirus   |
| <b>Target Names</b>      | 27  |
| <b>Protein Names</b>     | Recommended name: Baseplate structural protein Gp27 Alternative name(s):<br>Hub protein 27  |
| <b>Expression Region</b> | 1-391   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |