



Recombinant Enterobacteria phage T4 Baseplate structural protein Gp27 (27)

Product Code	CSB-MP325483EDZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P17172
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Purity	>85% (SDS-PAGE)
Sequence	MSMLQRPGYP NLSVKLFDSY DAWSNRNFVE LAATITTLTM RDSLYGRNEG MLQFYDSKNI HTKMDGNEII QISVANANDI NNVKTRIYGC KHFSVSVDSK GDNIIAIELG TIHSIENLKF GRPFFPDAGE SIKEMLGVIY QDRLLTPAI NAINAYVPDI PWTSTFENYL SYVREVALAV GSDKFVFWQ DIMGVNMMDY DMMINQEPYP MIVGEP SLIG QFIQELKYPL AYDFVWLT KS NPHKRDP MKN ATYAH SFLD SSIPMITTGK GENSIVV SRS GAYSEMTYRN GYEEAIRLQT MAQYDGYAKC STIGNFNLTP GVKIIFNSK NQFKTEFYVD EVIHLSNNN SVTHLYMFTN ATKLETIDPV KVKNEFKSDT TTESSSSNK Q
Source	Mammalian cell
Target Names	27
Protein Names	Recommended name: Baseplate structural protein Gp27 Alternative name(s): Hub protein 27
Expression Region	1-391
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.