



Recombinant *Saccharomyces cerevisiae* Malate dehydrogenase, mitochondrial (MDH1)

Product Code	CSB-YP325508SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P17505
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	YKV TVLGAGGGIG QPLSLLLKLN HKVTDLRLYD LKGAKGVATD LSHIPTNSVV KGFTPEEPDG LNNALKD TDM VLIPAGVPRK PGMTRDDLFA INASIVRDLA AATAESAPNA AILVISNPVN STVPIVAQVL KNKGVYNPKK LFGVTTLDSI RAARFISEVE NTDPTQERVN VIGGHSGITI IPLISQTNHK LMSDDKRHEL IHRIQFGGDE VVKAKNGAGS ATLSMAHAGA KFANAVLSGF KGERDVIEPS FVDSPLFKSE GIEFFASPVT LGPDGIEKIH PIGELSSEEE EMLQKCKETL KKNIEKGVNF VASK
Source	Yeast
Target Names	MDH1
Protein Names	Recommended name: Malate dehydrogenase, mitochondrial EC= 1.1.1.37
Expression Region	18-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.