



Recombinant Zoarces americanus Ice-structuring protein lambda OP-3

Product Code	CSB-BP325687MAQ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P19606
Product Type	Recombinant Protein
Immunogen Species	Zoarces americanus (Ocean pout) (Macrozoarces americanus)
Purity	>85% (SDS-PAGE)
Sequence	NQSVVATQ LIPINTALT L VMMTTRVIYP TGIPAEDIPR LVSMQVNQAV PMGTTLMPDM VKFYCLCAPK N
Source	Baculovirus
Protein Names	Recommended name: Ice-structuring protein lambda OP-3 Short name= ISP lambda OP-3 Alternative name(s): Antifreeze protein lambda OP-3
Expression Region	23-91
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.