



Recombinant *Saccharomyces cerevisiae* Isopentenyl-diphosphate Delta-isomerase (IDI1)

Product Code	CSB-BP325755SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P15496
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	TADNNSMPH GAVSSYAKLV QNQTPEIDLE EFPEIPLQQ RPNTRSSETS NDESGTCFS GHDEEQIKLM NENCIVLDWD DNAIGAGTKK VCHLMENIEK GLLHRAFSVF IFNEQGELLL QQRATEKITF PDLWTNTCCS HPLCIDDELG LKGKLDDKIK GAITAAVRKL DHELGIPEDE TKTRGKFHFL NRIHYMAPSN EPWGEHEIDY ILFYKINAKE NLTVPNPVNE VRDFKWVSPN DLKTMFADPS YKFTPWFKII CENYLFNWWL QLDDLSEVEN DRQIHRML
Source	Baculovirus
Target Names	IDI1
Protein Names	Recommended name: Isopentenyl-diphosphate Delta-isomerase EC= 5.3.3.2 Alternative name(s): Isopentenyl pyrophosphate isomerase Short name= IPP isomerase
Expression Region	2-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.